

GK Antibody, Biotin conjugated



PACO46433

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Key enzyme in the regulation of glycerol uptake and metabolism.

Gene ID:

GK

Uniprot

P32189

Synonyms:

Glycerol kinase (GK) (Glycerokinase) (EC 2.7.1.30) (ATP: glycerol 3-phosphotransferase), GK

Immunogen:

Recombinant Human Glycerol kinase protein (346-509AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A