

VACWR062 Antibody



PACO46438

Product Information

Size:

50ug

Reactivity:

Vaccinia virus

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Late protein which plays an essential role in virion assembly. Participates in the correct localization and aggregation of the viral core proteins.

Gene ID:

VACWR062

Uniprot

P21607

Synonyms:

Protein E6

Immunogen:

Recombinant Vaccinia virus Protein E6 protein (180-429AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A