

NRXN1 Antibody, Biotin conjugated



PACO46449

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Neuronal cell surface protein that may be involved in cell recognition and cell adhesion by forming intracellular junctions through binding to neuroligins. May play a role in formation or maintenance of synaptic junctions. May mediate intracellular signaling. May play a role in angiogenesis.

Gene ID:

NRXN1

Uniprot

P58400

Synonyms:

Neurexin-1-beta (Neurexin I-beta), NRXN1

Immunogen:

Recombinant Human Neurexin-1-β protein (184-363AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A