

Mup6 Antibody



PACO46450

Product Information

Size:

50ug

Reactivity:

Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IF

Recommended dilutions:

ELISA:1:2000-1:10000, IF:1:50-1:200

Protein Background:

Binds pheromones that are released from drying urine of males. These pheromones affect the sexual behavior of females.

Gene ID:

Mup6

Uniprot

P02762

Synonyms:

Major urinary protein 6 (MUP 6) (Alpha-2U-globulin) (Group 1, BS6) (allergen Mus m 1), Mup6


Immunogen:

Recombinant Mouse Major urinary protein 6 protein (19-180AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunofluorescent analysis of HeLa cells using PACO46450 at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).