

# KDM1A Antibody, HRP conjugated



PACO46471

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

KDM1A

**Uniprot**

O60341

**Synonyms:**

Lysine-specific histone demethylase 1A (EC 1) (BRAF35-HDAC complex protein BHC110) (Flavin-containing amine oxidase domain-containing protein 2), KDM1A, AOF2 KDM1 KIAA0601 LSD1

**Immunogen:**

Recombinant Human Lysine-specific histone demethylase 1A protein (100-150AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A