

KDM1A Antibody, FITC conjugated



PACO46472

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

KDM1A

Uniprot

O60341

Synonyms:

Lysine-specific histone demethylase 1A (EC 1) (BRAF35-HDAC complex protein BHC110) (Flavin-containing amine oxidase domain-containing protein 2), KDM1A, AOF2 KDM1 KIAA0601 LSD1

Immunogen:

Recombinant Human Lysine-specific histone demethylase 1A protein (100-150AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A