

PTPRS Antibody, Biotin conjugated



PACO46485

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Interacts with LAR-interacting protein LIP.1.

Gene ID:

PTPRS

Uniprot

Q13332

Synonyms:

Receptor-type tyrosine-protein phosphatase S (R-PTP-S) (EC 3.1.3.48) (Receptor-type tyrosine-protein phosphatase sigma) (R-PTP-sigma), PTPRS

Immunogen:

Recombinant Human Receptor-type tyrosine-protein phosphatase S protein (25-128AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A