

ME3 Antibody, Biotin conjugated



PACO46493

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

mitochondrion, cofactor binding, malate dehydrogenase (decarboxylating) (NADP+) activity, malic enzyme activity, aerobic respiration, malate metabolic process, oxidation-reduction process, oxygen metabolic process, pyruvate metabolic process.

Gene ID:

ME3

Uniprot

Q16798

Synonyms:

NADP-dependent malic enzyme, mitochondrial (NADP-ME) (EC 1.1.1.40) (Malic enzyme 3), ME3

Immunogen:

Recombinant Human NADP-dependent malic enzyme, mitochondrial protein (547-598AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A