

PDK4 Antibody



PACO46502

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Gene ID:

PDK4

Uniprot

Q16654

Synonyms:

[Pyruvate dehydrogenase (acetyl-transferring)] kinase isozyme 4, mitochondrial (EC 2.7.11.2) (Pyruvate dehydrogenase kinase isoform 4), PDK4, PDHK4


Immunogen:

Recombinant Human [Pyruvate dehydrogenase (acetyl-transferring)] kinase isozyme 4, mitochondrial protein (143-333AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human heart tissue using PACO46502 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human kidney tissue using PACO46502 at dilution of 1:100.