

# PDK4 Antibody, HRP conjugated



PACO46503

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

PDK4

**Uniprot**

Q16654

**Synonyms:**

[Pyruvate dehydrogenase (acetyl-transferring)] kinase isozyme 4, mitochondrial (EC 2.7.11.2) (Pyruvate dehydrogenase kinase isoform 4), PDK4, PDHK4

**Immunogen:**

Recombinant Human [Pyruvate dehydrogenase (acetyl-transferring)] kinase isozyme 4, mitochondrial protein (143-333AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A