

# STIL Antibody, HRP conjugated



PACO46515

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Immediate-early gene. Plays an important role in embryonic development as well as in cellular growth and proliferation; its long-term silencing affects cell survival and cell cycle distribution as well as decreases CDK1 activity correlated with reduced phosphorylation of CDK1. Plays a role as a positive regulator of the sonic hedgehog pathway, acting downstream of PTCH1.

**Gene ID:**

STIL

**Uniprot**

Q15468

**Synonyms:**

SCL-interrupting locus protein (TAL-1-interrupting locus protein), STIL, SIL

**Immunogen:**

Recombinant Human SCL-interrupting locus protein (388-637AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A