

# KYNU Antibody, FITC conjugated



PACO46548

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Catalyzes the cleavage of L-kynurenine (L-Kyn) and L-3-hydroxykynurenine (L-3OHKyn) into anthranilic acid, (AA) and 3-hydroxyanthranilic acid, (3-OHAA), respectively. Has a preference for the L-3-hydroxy form. Also has cysteine-conjugate-beta-lyase activity.

**Gene ID:**

KYNU

**Uniprot**

Q16719

**Synonyms:**

Kynureninase (EC 3.7.1.3) (L-kynurenine hydrolase), KYNU

**Immunogen:**

Recombinant Human Kynureninase protein (99-214AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A