

JMJD1C Antibody, FITC conjugated



PACO46568

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Probable histone demethylase that specifically demethylates 'Lys-9' of histone H3, thereby playing a central role in histone code. Demethylation of Lys residue generates formaldehyde and succinate. May be involved in hormone-dependent transcriptional activation, by participating in recruitment to androgen-receptor target genes.

Gene ID:

JMJD1C

Uniprot

Q15652

Synonyms:

Probable JmjC domain-containing histone demethylation protein 2C (EC 1.14.11) (Jumonji domain-containing protein 1C) (Thyroid receptor-interacting protein 8) (TR-interacting protein 8) (TRIP-8), JMJD1C, JHDM2C KIAA1380 TRIP8

Immunogen:

Recombinant Human Probable JmjC domain-containing histone demethylation protein 2C protein (539-648AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A