

ACSF3 Antibody, HRP conjugated



PACO46579

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Catalyzes the initial reaction in intramitochondrial fatty acid, synthesis, by activating malonate and methylmalonate, but not acetate, into their respective CoA thioester. May have some preference toward very-long-chain substrates.

Gene ID:

ACSF3

Uniprot

Q4G176

Synonyms:

Acyl-CoA synthetase family member 3, mitochondrial (EC 6.2.1), ACSF3

Immunogen:

Recombinant Human Acyl-CoA synthetase family member 3, mitochondrial protein (365-484AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A