

TSBP1 Antibody, HRP conjugated



PACO46591

Product Information

Size:	Protein Background:
50ug	nucleus.
Reactivity:	Gene ID:
Human	TSBP1
Source:	Uniprot
Rabbit	Q5SRN2
Isotype:	Synonyms:
IgG	Testis-expressed basic protein 1(Uncharacterized protein C6orf10), TSBP1, C6orf10
Applications:	Immunogen:
ELISA	Recombinant Human Testis-expressed basic protein 1 protein (239-476AA).
Recommended dilutions:	Storage:
	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A