

# TSBP1 Antibody, Biotin conjugated



PACO46593

---

## Product Information

<b>Size:</b>	<b>Protein Background:</b>
50ug	nucleus.
<b>Reactivity:</b>	<b>Gene ID:</b>
Human	TSBP1
<b>Source:</b>	<b>Uniprot</b>
Rabbit	Q5SRN2
<b>Isotype:</b>	<b>Synonyms:</b>
IgG	Testis-expressed basic protein 1(Uncharacterized protein C6orf10), TSBP1, C6orf10
<b>Applications:</b>	<b>Immunogen:</b>
ELISA	Recombinant Human Testis-expressed basic protein 1 protein (239-476AA).
<b>Recommended dilutions:</b>	<b>Storage:</b>
	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A