

RUNX1T1 Antibody, FITC conjugated



PACO46608

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Transcriptional corepressor which facilitates transcriptional repression via its association with DNA-binding transcription factors and recruitment of other corepressors and histone-modifying enzymes. Can repress the expression of MMP7 in a ZBTB33-dependent manner. Can repress transactivation mediated by TCF12. Acts as a negative regulator of adipogenesis.

Gene ID:

RUNX1T1

Uniprot

Q06455

Synonyms:

Protein CBFA2T1 (Cyclin-D-related protein) (Eight twenty one protein) (Protein ETO) (Protein MTG8) (Zinc finger MYND domain-containing protein 2), RUNX1T1, AML1T1 CBFA2T1 CDR ETO MTG8 ZMYND2

Immunogen:

Recombinant Human Protein CBFA2T1 protein (416-525AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A