

RABGAP1L Antibody, Biotin conjugated



PACO46625

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

RABGAP1L

Uniprot

B7ZAP0

Synonyms:

Rab GTPase-activating protein 1-like, isoform 10, RABGAP1L, HHL KIAA0471

Immunogen:

Recombinant Human Rab GTPase-activating protein 1-like, isoform 10 protein (11-262AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A