

# TENT5C Antibody, HRP conjugated



PACO46635

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Seems to enhance replication of some viruses, including yellow fever virus, in response to type I interferon.

**Gene ID:**

TENT5C

**Uniprot**

Q5VWP2

**Synonyms:**

Putative nucleotidyltransferase FAM46C (EC 2.7.7), FAM46C

**Immunogen:**

Recombinant Human Terminal nucleotidyltransferase 5C protein (1-391AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A