

CDC73 Antibody



PACO46678

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IF

Recommended dilutions:

ELISA:1:2000-1:10000, IF:1:50-1:200

Protein Background:

Gene ID:

CDC73

Uniprot

Q6P1J9

Synonyms:

Parafibromin (Cell division cycle protein 73 homolog) (Hyperparathyroidism 2 protein), CDC73, C1orf28 HRPT2


Immunogen:

Recombinant Human Parafibromin protein (74-346AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunofluorescent analysis of HeLa cells using PACO46678 at dilution of 1: 100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).