

CDC73 Antibody, HRP conjugated



PACO46679

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

CDC73

Uniprot

Q6P1J9

Synonyms:

Parafibromin (Cell division cycle protein 73 homolog) (Hyperparathyroidism 2 protein), CDC73, C1orf28 HRPT2

Immunogen:

Recombinant Human Parafibromin protein (74-346AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A