

NBEAL2 Antibody, FITC conjugated



PACO46688

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Probably involved in thrombopoiesis. Plays a role in the development or secretion of alpha-granules, that contain several growth factors important for platelet biogenesis.

Gene ID:

NBEAL2

Uniprot

Q6ZNJ1

Synonyms:

Neurobeachin-like protein 2, NBEAL2, KIAA0540

Immunogen:

Recombinant Human Neurobeachin-like protein 2 protein (1360-1574AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A