

FGD5 Antibody, FITC conjugated



PACO46700

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Activates CDC42, a member of the Ras-like family of Rho- and Rac proteins, by exchanging bound GDP for free GTP. Mediates VEGF-induced CDC42 activation. May regulate proangiogenic action of VEGF in vascular endothelial cells, including network formation, directional movement and proliferation. May play a role in regulating the actin cytoskeleton and cell shape.

Gene ID:

FGD5

Uniprot

Q6ZNL6

Synonyms:

FYVE, RhoGEF and PH domain-containing protein 5 (Zinc finger FYVE domain-containing protein 23), FGD5, ZFYVE23

Immunogen:

Recombinant Human FYVE, RhoGEF and PH domain-containing protein 5 protein (269-539AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A