

DET1 Antibody, FITC conjugated



PACO46704

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Component of the E3 ubiquitin ligase DCX DET1-COP1 complex, which is required for ubiquitination and subsequent degradation of target proteins. The complex is involved in JUN ubiquitination and degradation.

Gene ID:

DET1

Uniprot

Q7L5Y6

Synonyms:

DET1 homolog (De-etiolated-1 homolog), DET1

Immunogen:

Recombinant Human DET1 homolog protein (5-190AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A