

PRMT9 Antibody



PACO46706

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Gene ID:

PRMT9

Uniprot

Q6P2P2

Synonyms:

Protein arginine N-methyltransferase 9 (Protein arginine N-methyltransferase 10) (EC 2.1.1.320), PRMT9, PRMT10


Immunogen:

Recombinant Human Putative protein arginine N-methyltransferase 9 protein (181-453AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human kidney tissue using PACO46706 at dilution of 1:100.