

KANSL1 Antibody, FITC conjugated



PACO46716

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

As part of the NSL complex it is involved in acetylation of nucleosomal histone H4 on several lysine residues and therefore may be involved in the regulation of transcription.

Gene ID:

KANSL1

Uniprot

Q7Z3B3

Synonyms:

KAT8 regulatory NSL complex subunit 1 (MLL1/MLL complex subunit KANSL1) (MSL1 homolog 1) (hMSL1v1) (NSL complex protein NSL1) (Non-specific lethal 1 homolog), KANSL1, CENP-36 KIAA1267 MSL1V1 NSL1

Immunogen:

Recombinant Human KAT8 regulatory NSL complex subunit 1 protein (815-1086AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A