

MORC1 Antibody, HRP conjugated



PACO46751

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Required for spermatogenesis.

Gene ID:

MORC1

Uniprot

Q86VD1

Synonyms:

MORC family CW-type zinc finger protein 1 (Cancer/testis antigen 33) (CT33), MORC1, MORC

Immunogen:

Recombinant Human MORC family CW-type zinc finger protein 1 protein (683-984AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A