

ALDH16A1 Antibody, FITC conjugated



PACO46776

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

extracellular exosome, membrane, 4-trimethylammoniobutyraldehyde dehydrogenase activity, aminobutyraldehyde dehydrogenase activity.

Gene ID:

ALDH16A1

Uniprot

Q8IZ83

Synonyms:

Aldehyde dehydrogenase family 16 member A1, ALDH16A1

Immunogen:

Recombinant Human Aldehyde dehydrogenase family 16 member A1 protein (1-220AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A