

FNIP1 Antibody, FITC conjugated



PACO46812

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

May be involved in energy and/or nutrient sensing through the AMPK and mTOR signaling pathways. May regulate phosphorylation of RPS6KB1.

Gene ID:

FNIP1

Uniprot

Q8TF40

Synonyms:

Folliculin-interacting protein 1, FNIP1, KIAA1961

Immunogen:

Recombinant Human Folliculin-interacting protein 1 protein (221-508AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A