

RTF1 Antibody, HRP conjugated



PACO46819

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

RTF1

Uniprot

Q92541

Synonyms:

RNA polymerase-associated protein RTF1 homolog, RTF1, KIAA0252

Immunogen:

Recombinant Human RNA polymerase-associated protein RTF1 homolog protein (620-670AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A