

PDLIM5 Antibody, HRP conjugated



PACO46827

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

May play an important role in the heart development by scaffolding PKC to the Z-disk region. May play a role in the regulation of cardiomyocyte expansion. Overexpression promotes the development of heart hypertrophy. Contributes to the regulation of dendritic spine morphogenesis in neurons. May restrain postsynaptic growth of excitatory synapses.

Gene ID:

PDLIM5

Uniprot

Q96HC4

Synonyms:

PDZ and LIM domain protein 5 (Enigma homolog) (Enigma-like PDZ and LIM domains protein), PDLIM5, ENH

Immunogen:

Recombinant Human PDZ and LIM domain protein 5 protein (268-404AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A