

Nlrp3 Antibody, HRP conjugated



PACO46847

Product Information

Size:	Protein Background:
50ug	
Reactivity:	Gene ID:
Mouse	Nlrp3
Source:	Uniprot
Rabbit	Q8R4B8
Isotype:	Synonyms:
IgG	NACHT, LRR and PYD domains-containing protein 3 (Cold autoinflammatory syndrome 1 protein homolog) (Cryopyrin) (Mast cell maturation-associated-inducible protein 1) (PYRIN-containing APAF1-like protein 1), Nlrp3, Cias1 Mmig1 Nalp3 Pypaf1
Applications:	Immunogen:
ELISA	Recombinant Mouse NACHT, LRR and PYD domains-containing protein 3 protein (1-153AA).
Recommended dilutions:	Storage:
	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A