

# P4HB Antibody, FITC conjugated



PACO46852

---

## Product Information

**Size:**

50ug

**Reactivity:**

Cricetulus griseus

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

P4HB

**Uniprot**

Q8R4U2

**Synonyms:**

Protein disulfide-isomerase (PDI) (EC 5.3.4.1) (Prolyl 4-hydroxylase subunit beta) (p58), P4HB, PDIA1

**Immunogen:**

Recombinant Cricetulus griseus Protein disulfide-isomerase protein (20-509AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A