

P4HB Antibody, Biotin conjugated



PACO46853

Product Information

Size:

50ug

Reactivity:

Cricetulus griseus

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

P4HB

Uniprot

Q8R4U2

Synonyms:

Protein disulfide-isomerase (PDI) (EC 5.3.4.1) (Prolyl 4-hydroxylase subunit beta) (p58), P4HB, PDIA1

Immunogen:

Recombinant Cricetulus griseus Protein disulfide-isomerase protein (20-509AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A