

# RNF125 Antibody, HRP conjugated



PACO46863

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

E3 ubiquitin-protein ligase that acts as a positive regulator of T-cell activation. E3 ligase proteins mediate ubiquitination and subsequent proteasomal degradation of target proteins.

**Gene ID:**

RNF125

**Uniprot**

Q96EQ8

**Synonyms:**

E3 ubiquitin-protein ligase RNF125 (EC 2.3.2.27) (RING finger protein 125) (T-cell RING activation protein 1) (TRAC-1), RNF125

**Immunogen:**

Recombinant Human E3 ubiquitin-protein ligase RNF125 protein (143-231AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A