

RNF125 Antibody, Biotin conjugated



PACO46865

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

E3 ubiquitin-protein ligase that acts as a positive regulator of T-cell activation. E3 ligase proteins mediate ubiquitination and subsequent proteasomal degradation of target proteins.

Gene ID:

RNF125

Uniprot

Q96EQ8

Synonyms:

E3 ubiquitin-protein ligase RNF125 (EC 2.3.2.27) (RING finger protein 125) (T-cell RING activation protein 1) (TRAC-1), RNF125

Immunogen:

Recombinant Human E3 ubiquitin-protein ligase RNF125 protein (143-231AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A