

ACTL9 Antibody, HRP conjugated



PACO46895

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

ACTL9

Uniprot

Q8TC94

Synonyms:

Actin-like protein 9, ACTL9, ACTL7C

Immunogen:

Recombinant Human Actin-like protein 9 protein (1-160AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A