

AVL9 Antibody



PACO46914

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC, IF

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200,
IF:1:50-1:200

Protein Background:

Functions in cell migration.

Gene ID:

AVL9

Uniprot

Q8NBF6

Synonyms:

Late secretory pathway protein AVL9 homolog, AVL9, KIAA0241


Immunogen:

Recombinant Human Late secretory pathway protein AVL9 homolog protein (162-432AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human pancreatic tissue using PACO46914 at dilution of 1:100.



Immunofluorescent analysis of HepG2 cells using PACO46914 at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).