

THRSP Antibody, Biotin conjugated



PACO46933

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Plays a role in the regulation of lipogenesis, especially in lactating mammary gland. Important for the biosynthesis of triglycerides with medium-length fatty acid, chains. May modulate lipogenesis by interacting with MID1IP1 and preventing its interaction with ACACA. May function as transcriptional coactivator. May modulate the transcription factor activity of THRβ.

Gene ID:

THRSP

Uniprot

Q92748

Synonyms:

Thyroid hormone-inducible hepatic protein (Spot 14 protein) (S14) (SPOT14), THRSP

Immunogen:

Recombinant Human Thyroid hormone-inducible hepatic protein (41-144AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A