

GATA6 Antibody



PACO46934

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Transcriptional activator that regulates SEMA3C and PLXNA2. Involved in gene regulation specifically in the gastric epithelium. May regulate genes that protect epithelial cells from bacterial infection. Involved in bone morphogenetic protein (BMP)-mediated cardiac-specific gene expression. Binds to BMP response element (BMPRE) DNA sequences within cardiac activating regions.

Gene ID:

GATA6

Uniprot

Q92908

Synonyms:

Transcription factor GATA-6 (GATA-binding factor 6), GATA6


Immunogen:

Recombinant Human Transcription factor GATA-6 protein (496-595AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human brain tissue using PACO46934 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human melanoma using PACO46934 at dilution of 1:100.