

GATA6 Antibody, HRP conjugated



PACO46935

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Transcriptional activator that regulates SEMA3C and PLXNA2. Involved in gene regulation specifically in the gastric epithelium. May regulate genes that protect epithelial cells from bacterial infection. Involved in bone morphogenetic protein (BMP)-mediated cardiac-specific gene expression. Binds to BMP response element (BMPRE) DNA sequences within cardiac activating regions.

Gene ID:

GATA6

Uniprot

Q92908

Synonyms:

Transcription factor GATA-6 (GATA-binding factor 6), GATA6

Immunogen:

Recombinant Human Transcription factor GATA-6 protein (496-595AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A