

# DNAAF3 Antibody



PACO46962

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC, IF

**Recommended dilutions:**

ELISA:1:2000-1:10000, IHC:1:20-1:200,  
IF:1:50-1:200

**Protein Background:**

Required for the assembly of axonemal inner and outer dynein arms. Involved in preassembly of dyneins into complexes before their transport into cilia.

**Gene ID:**

DNAAF3

**Uniprot**

Q8N9W5

**Synonyms:**

Dynein assembly factor 3, axonemal, DNAAF3, C19orf51

**Immunogen:**


Recombinant Human Dynein assembly factor 3, axonemal protein (1-209AA).

**Storage:**


Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---



Immunohistochemistry of paraffin-embedded human appendix tissue using PACO46962 at dilution of 1:100.



Immunofluorescent analysis of HepG2 cells using PACO46962 at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).