

PRMT6 Antibody



PACO47006

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC, IF

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200,
IF:1:50-1:200

Protein Background:**Gene ID:**

PRMT6

Uniprot

Q96LA8

Synonyms:

Protein arginine N-methyltransferase 6 (EC 2.1.1.319) (Heterogeneous nuclear ribonucleoprotein methyltransferase-like protein 6) (Histone-arginine N-methyltransferase PRMT6), PRMT6, HRMT1L6


Immunogen:

Recombinant Human Protein arginine N-methyltransferase 6 protein (119-375AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4


Product Images



Immunohistochemistry of paraffin-embedded human testis tissue using PACO47006 at dilution of 1:100.



Immunofluorescent analysis of HeLa cells using PACO47006 at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).



Immunohistochemistry of paraffin-embedded human kidney tissue using PACO47006 at dilution of 1:100.