

SCPEP1 Antibody, Biotin conjugated



PACO47065

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

May be involved in vascular wall and kidney homeostasis.

Gene ID:

SCPEP1

Uniprot

Q9HB40

Synonyms:

Retinoid-inducible serine carboxypeptidase (EC 3.4.16) (Serine carboxypeptidase 1), SCPEP1, RISC SCP1

Immunogen:

Recombinant Human Retinoid-inducible serine carboxypeptidase protein (201-452AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A