

SMOC1 Antibody, Biotin conjugated



PACO47105

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Plays essential roles in both eye and limb development. Probable regulator of osteoblast differentiation.

Gene ID:

SMOC1

Uniprot

Q9H4F8

Synonyms:

SPARC-related modular calcium-binding protein 1 (Secreted modular calcium-binding protein 1) (SMOC-1), SMOC1

Immunogen:

Recombinant Human SPARC-related modular calcium-binding protein 1 protein (277-382AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A