

UVRAG Antibody, FITC conjugated



PACO47156

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

UVRAG

Uniprot

Q9P2Y5

Synonyms:

UV radiation resistance-associated gene protein (p63), UVRAG

Immunogen:

Recombinant Human UV radiation resistance-associated gene protein (375-425AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A