

# SLC25A10 Antibody, Biotin conjugated



PACO47165

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Involved in translocation of malonate, malate and succinate in exchange for phosphate, sulfate, sulfite or thiosulfate across mitochondrial inner membrane.

**Gene ID:**

SLC25A10

**Uniprot**

Q9UBX3

**Synonyms:**

Mitochondrial dicarboxylate carrier (Solute carrier family 25 member 10), SLC25A10, DIC

**Immunogen:**

Recombinant Human Mitochondrial dicarboxylate carrier protein (124-162AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A