

BTNL2 Antibody, FITC conjugated



PACO47172

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Negative regulator of T-cell proliferation.

Gene ID:

BTNL2

Uniprot

Q9UIR0

Synonyms:

Butyrophilin-like protein 2 (BTL-II), BTNL2

Immunogen:

Recombinant Human Butyrophilin-like protein 2 protein (26-202AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A