

# BTNL2 Antibody, Biotin conjugated



PACO47173

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Negative regulator of T-cell proliferation.

**Gene ID:**

BTNL2

**Uniprot**

Q9UIR0

**Synonyms:**

Butyrophilin-like protein 2 (BTL-II), BTNL2

**Immunogen:**

Recombinant Human Butyrophilin-like protein 2 protein (26-202AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A