

EHD1 Antibody, HRP conjugated



PACO47183

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

EHD1

Uniprot

Q9H4M9

Synonyms:

EH domain-containing protein 1 (PAST homolog 1) (hPAST1) (Testilin), EHD1, PAST PAST1

Immunogen:

Recombinant Human EH domain-containing protein 1 protein (401-498AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A